Baked Mozzarella Bites

Cooking Light Magazine - August 2011

Servings: 4

1/3 cup panko (Japanese breadcrumbs)
3 sticks (1 ounce ea) part-skim mozzarella string cheese
3 tablespoons egg substitute
cooking spray
1/4 cup lower-sodium marinara sauce

Preheat the oven to 425 degrees.

Heat a medium skillet over medium heat.

Add the panko to the pan. Cook for 2 minutes or until toasted, stirring frequently. Remove from the heat and place the panko in a shallow dish.

Cut the mozzarella sticks into 1-inch pieces.

Working with one piece at a time, dip the cheese into the egg substitute and dredge in the panko.

Place the cheese on a baking sheet coated with cooking spray.

Bake for 3 minutes or until the cheese is softened and thoroughly heated.

Pour the marinara sauce into a microwave-safe bowl. Microwave on HIGH for 1 minute or until thoroughly heated, stirring after 30 seconds.

Serve with the mozzarella pieces.

Yield: 12 mozzarella bites

Per Serving (excluding unknown items): 18 Calories; 1g Fat (63.3% calories from fat); 1g Protein; trace Carbohydrate; 0g Dietary Fiber; trace Cholesterol; 22mg Sodium. Exchanges: 0 Lean Meat; 0 Fat.